


SURG-I-BAND®

COLOR CODING



Original SCANLAN® Surg-I-Band® (1/4") and narrow SCANLAN® Surg-I-Band® (1/8") color coding, available in 14 bright colors, offer a multitude of color coding combinations to identify instruments by surgeon, set, specialty and suite. In addition, Surg-I-Band® color coding is long lasting – providing greater value by saving time and labor. This high quality product has been independently tested and confirmed to be completely gas and steam sterilizable. Not made with natural rubber latex.

1/8" Narrow Cat. #	1/4" Original Cat. #	Description	1/8" Narrow Cat. #	1/4" Original Cat. #	Description
1001-630	1001-64	Red	1001-637	1001-71	Black
1001-631	1001-65	Yellow	1001-638	1001-94	Glo-Orange
1001-632	1001-66	Orange	1001-639	1001-95	Purple
1001-633	1001-67	Blue	1001-640	1001-96	Lime Green
1001-634	1001-68	Green	1001-641	1001-97	Grey
1001-635	1001-69	Brown	1001-642	1001-98	Pink
1001-636	1001-70	White	1001-643	1001-99	Light Blue

SURG-I-BAND® Dispenser



- **Easy to use** – SCANLAN® Surg-I-Band® color coding spools snap on securely and can be removed quickly
- **Precise cut** – Dispenser eliminates the need for awkward scissors or blades when using SCANLAN® Surg-I-Band® color coding
- **Protective cover** - guards against damage to blade and color coding
- **Compact** – provides a small, lightweight design
- **Reusable** – ensuring greater economy

Cat. #	Description
1001-60	Surg-I-Band® Dispenser (Does not include Surg-I-Band® spools)

SURG-I-BAND® Remover



U.S. Patent Number 5,864,953

- **Saves time** – allows quick removal of SCANLAN® Surg-I-Band® color coding
- **Safe** – eliminates use of razors or other sharp instruments that may damage the surface of instruments or cause injury to personnel
- **Reusable** – removes hundreds of bands

Cat. #	Description
1001-49	Surg-I-Band® Remover

SURG-I-BAND® Carousel



- **Efficient** – provides fast, easy application of SCANLAN® Surg-I-Band® color coding
- **Excellent organization** – holds 12 SCANLAN® Surg-I-Band® dispensers
- **Rotating top** – provides easy access to each color
- **Storage** – conveniently holds Scanlan® Tip Guard™ instrument protectors
- **Easy to use** – SCANLAN® Surg-I-Band® spools snap on securely and can be removed quickly

Cat. #	Description
1001-59	Carousel Only
1001-61	Carousel with 12 Dispensers (Does not include Surg-I-Band® spools or Tip-Guards™)

