







Portfolio Overview

Original-Perfusor® Lines

	Description	Tube material	Length (cm)	Ø (mm)	Priming vol. (ml)	Units per box	Code no. (REF)
 <ul style="list-style-type: none"> Luer-Lock patient connector DEHP-free 		PVC	25	1.5 x 2.7	0.5	100	0066088K
		PVC	50	1.5 x 2.7	1.3	100	8255172
		PVC	100	1.5 x 2.7	2.0	100	0009483H
		PVC	150	1.5 x 2.7	2.9	100	8722960
		PVC	200	1.5 x 2.7	3.8	100	8722862
		PVC	250	1.5 x 2.7	4.8	100	8255490
		PVC	300	1.5 x 2.7	5.7	100	8255253
		PVC	150	1.5 x 2.7	3.3	100	8722820*
		PVC	200	1.5 x 2.7	4.8	100	8723001**
		PVC	75	0.9 x 1.9	0.6	100	8722870N
		PVC	150	0.9 x 1.9	1.1	100	8255504N
		PVC	200	0.9 x 1.9	1.4	100	8745919N
 <ul style="list-style-type: none"> Anti-Free-Flow Protection Clamp Luer-Lock patient connector DEHP-free 		PVC	150	0.9 x 1.9	1.1	100	8722940
		PVC	200	0.9 x 1.9	1.4	100	8722941

* with Safsite® Valve

** with Sterifix® 0.2 µ

	Description	Tube material	Length (cm)	Ø (mm)	Priming vol. (ml)	Units per box	Code no. (REF)
 <ul style="list-style-type: none"> Luer-Lock patient connector PVC-free 		PE	10	1.0 x 2.0	0.2	100	0066086H
		PE	25	1.0 x 2.0	0.3	100	0066087J
		PE	50	1.0 x 2.0	0.5	100	8255059
		PE	100	1.0 x 2.0	0.9	100	8255067
		PE	150	1.0 x 2.0	1.3	100	8722935
		PE	200	1.0 x 2.0	1.7	100	8723060
		PE	250	1.0 x 2.0	2.1	100	8272565
		PE	300	1.0 x 2.0	2.5	100	8250146
 <ul style="list-style-type: none"> Light protection up to 520 nm Luer-Lock patient connector PVC-free 		PE	150	1.0 x 2.0	1.3	100	8723017
		PE	200	1.0 x 2.0	1.7	100	8723018
 <ul style="list-style-type: none"> ENFit™ female; ENFit™ male 		PE	150	2.5 x 4.3	7.9	25	87229910
 <ul style="list-style-type: none"> NRFit® patient connector Translucent yellow tube 		PUR	200	1.3 x 2.4	2.8	25	8722865