

## Cryoalfa<sup>®</sup> SUPER and SUPER-CONTACT

The most popular and simplest device for daily use by medical professionals and beauticians.

**SUPER** The freezing tip is fixed and not interchangeable.

**SUPER CONTACT** The freezing tip is not interchangeable.

There are 4 devices with different applicators to choose from.  
Use of cartridges with valve.

**CA-S Super - Liquid Freezing** glass tip Ø 1mm

**CA-S-C Super - Contact Freezing**

**CA-S-C3** gold-plated contact tip Ø 3mm

**CA-S-C5** gold-plated contact surface Ø 5mm

**CA-S-C7** gold-plated contact surface Ø 7mm

