


# Sharp Needles



Hub Features	PRE	PRC	HPC	LDS	Hub Benefits
Hard Pack Capsule	X	X	X	X	Preventing Needle Injury
Hard Polymer		X	X	X	No Flexing
Threading			X	X	No Pop-Off Risk
Low Dead Space				X	No Product Loss
Needle Up to 14% thinner than 33G				X	Minimal Patient Trauma
Recommended Use	BoNT + Meso	Filler	Filler	BoNT	

## THE INViSIBLE NEEDLE™.


**THE INViSIBLE NEEDLE™** is the first needle of its kind. Technology advancements have made it possible to make this needle 33% thinner than a 30G needle (used for insulin syringes) and 14% thinner than the TSK 33G needle. This extremely small diameter of THE INViSIBLE NEEDLE™ results in a mere painless experience for patients and is ideal for BoNT injections.



**NEW**

## Low dead space hub.

**The Low Dead Space needle hub** is designed to reduce the effective space left behind in the needle hub as close to zero. This new patented LDS hub design results in savings of up to 0.08ml which leads to significant cost savings when injecting high cost drugs, such as botulinum toxin.



### BoNT Needles / LDS Hub

Code	Description	Inch	UTW*	Box
LDS-02009	THE INViSIBLE NEEDLE x 9mm**	3/8"		100
LDS-33013	Low Dead Space Hub 33G x 13mm**	1/2"		100
LDS-33009	Low Dead Space Hub 33G x 9mm**	3/8"		100

### BoNT & Mesotherapy Needles / PRE Hub

Code	Description	Inch	UTW*	Box
PRE-33013	Regular Hub 33G x 13mm	1/2"		100
PRE-33004	Regular Hub 33G x 4mm	3/16"		100
PRE-32013	Regular Hub 32G x 13mm	1/2"		100
PRE-32009	Regular Hub 32G x 9mm	3/8"		100
PRE-32006	Regular Hub 32G x 6mm	1/4"		100
PRE-32004	Regular Hub 32G x 4mm	3/16"		100
PRE-31013	Regular Hub 31G x 13mm	1/2"		100
PRE-30006I	Regular Hub 30G x 6mm	1/4"	✓	100
PRE-30004I	Regular Hub 30G x 4mm	3/16"	✓	100

### Filler Needles / PRC Hub

Code	Description	Inch	UTW*	Box
PRC-30038I	Control Hub 30G x 38mm	1 1/2"	✓	100
PRC-30025I	Control Hub 30G x 25mm	1"	✓	100
PRC-30013I	Control Hub 30G x 13mm	1/2"	✓	100
PRC-27038I	Control Hub 27G x 38mm	1 1/2"	✓	100
PRC-27019I	Control Hub 27G x 19mm	3/4"	✓	100
PRC-27013I	Control Hub 27G x 13mm	1/2"	✓	100
PRC-25013I	Control Hub 25G x 13mm	1/2"	✓	100
PRC-22050I	Control Hub 22G x 50mm	2"	✓	100

### Filler Needles / HPC Hub

Code	Description	Inch	UTW*	Box
HPC-30013I	Advanced Hub 30G x 13mm	1/2"	✓	320
HPC-27013I	Advanced Hub 27G x 13mm	1/2"	✓	320
HPC-23025I	Advanced Hub 23G x 25mm	1"	✓	320

\* UTW = ULTRA Thin Wall \*\* For exclusive use with Soft-Ject Luer Lock syringe

All other sizes on request