

# Adapters and holders

## BD Vacutainer® single use holder, BD luer adapter and adapters with pre-attached holders

**1** BD Vacutainer® Luer adapter (Cat. No. 367300) is a sterile device with a multi-sampling valve and is designed for use with a catheter to collect blood with BD Vacutainer® blood collection tubes. The device is ideal for use with any female luer fitting.


**2** BD Vacutainer® single-use holders (Cat. No. 364815) are compatible with all BD Vacutainer® tubes and needles, including BD Eclipse™ safety needles, BD Safety-Lok™ blood collection sets and BD Push Button blood collection sets. BD Vacutainer® single-use holders are also compatible with BD BACTEC™ blood culture bottles.

**3** The BD Vacutainer® Luer-Lok™ Access Device (Cat. No. 364902) is a preassembled multisample BD Luer-Lok™ and holder which is compatible with female luer connections or IV ports designed for luer access and has a blue colour-coded connection to provide easy differentiation from other holder based-products.





**4** The BD Vacutainer® Blood Transfer Device (Cat. No. 364810) is a pre-assembled and easy-to-use device designed with safety in mind. It is used for needleless specimen transfer from a syringe to an evacuated tube or blood culture bottle and has a red colour-coded connection to provide easy differentiation from other holder based products.

## BD Vacutainer® Luer adapter




Cat No.	Description	Colour code	Unit of sales
367300	Luer extension for taking blood samples from catheters and perfusion sets		1000 (10x100)

## BD Vacutainer® Luer adapters with pre-attached holders

These single use products are ready-to-use, sterile, individually blister packaged holders, with the Luer adapter ready fitted.

Cat No.	Description	Colour code	Unit of sales
364902	BD Vacutainer® Luer-Lok™ Access Device ("male Luer")		200 (2x100)
364810	BD Vacutainer® Blood Transfer Device ("female Luer")		200 (2x100)

## BD Vacutainer® holders

Cat No.	Description	Colour code	Unit of sales
364815	BD Vacutainer® single use plastic holder for tubes with 13mm and 16mm diameter and for BD Bactec™ blood culture bottles		1000 (4x250)
368872	BD Pronto™ Quick Release holder for tubes with 13mm and 16mm diameter and for BD Bactec™ blood culture bottles	 various	100 (5x20)
364879	BD Vacutainer® Multiple use plastic holder for tubes with 13mm and 16mm diameter and for BD Bactec™ blood culture bottles		1000 (4x250)