

## BD™ Quincke spinal needle

- Key/slot arrangement of stylet and cannula hubs facilitates aligned needle and stylet bevel orientation
- Translucent window hub (featuring contact clarity) that helps visualization of cerebrospinal fluid
- Fitted stylet designed to reduce tissue coring
- Color-coded
- Sterile, single use



Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box/ Shipping Case	Reference
27 G	3 1/2"	90	Grey	25	405259
26 G	3 1/2"	90	Brown	25	405258
25 G	4 11/16"	119	Blue	10	405234
25 G	3 1/2"	90	Orange	25	405257
25 G	3"	75	Orange	25	405246
25 G	2"	51	Blue	25	405245
25 G	1"	25	Blue	25	405243
23 G	3 1/2"	90	Turquoise	25	405240
22 G	7"	178	Black	10	405149
22 G	5"	127	Black	10	405148
22 G	3 1/2"	90	Black	25	405256
22 G	3"	75	Black	25	405255
22 G	2 1/2"	63	Black	25	405244
22 G	1 1/2"	38	Black	25	405254
20 G	6"	152	Yellow	10	405211
20 G	3 1/2"	90	Yellow	25	405253
20 G	3"	75	Yellow	25	405252
20 G	1 1/2"	38	Yellow	25	405251
19 G	3 1/2"	90	Cream	25	405250
19 G	3"	75	Cream	25	405249
18 G	3 1/2"	90	Pink	25	405248
18 G	3"	75	Pink	25	405247

## BD™ Quincke spinal needle

With Introducer 22 G 1 1/4"

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
27 G	3 1/2"	90	Grey	25	405069

With Introducer 20 G 1 1/4"

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
26 G	3 1/2"	90	Brown	25	405065
25 G	3 1/2"	90	Orange	25	405084

BD™ Quincke Spinal Needle introducers

Gauge	Length (inch)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference
20 G	1 1/4"	32	Yellow	25	405260
22 G	1 1/4"	32	Black	25	405261