



IVALON


Nasal Packing

IVALON® Anatomical Nasal Packings

IVALON® Anatomical Nasal Packings are constructed of PVAc with the 4000Plus non-stick surface.

- Integral drawstring facilitates removal

Product Code	Product Description	Per Box
Q770409	Anatomical Nasal Packing, 4000Plus, Sterile 8.0 cm long x 1.5 cm wide x 3.0 cm high, 1 /Pkg; 10 Pkgs/Box	10
Q770411	Anatomical ThinPack™ Nasal Packing, 4000Plus, Sterile 8.0 cm long x 1.0 cm wide x 3.0 cm high, 1 /Pkg; 10 Pkgs/Box	10
Q771319	Anatomical Nasal Packing, No String, 4000Plus, Sterile 8.0 cm long x 1.5 cm wide x 3.0 cm high, 1 /Pkg; 20 Pkgs/Box	20




Anatomical Nasal Packings with Airway Tube

IVALON® Anatomical Nasal Packings with Airway Tube are constructed of PVAc with the 4000Plus non-stick surface.

- Integral drawstring facilitates removal

Product Code	Product Description	Per Box
Q770410	Anatomical Nasal Packing with Airway Tube, 4000Plus, Sterile 8.0 cm long x 1.5 cm wide x 3.0 cm high, 1 /Pkg; 10 Pkgs/Box	10




IVALON® Nasal Packings

IVALON® Nasal Packings are constructed of PVAc with the 4000Plus non-stick surface.

- Integral drawstring facilitates removal

Product Code	Product Description	Per Box
Q770401	Nasal Packing, 4000Plus, sterile 4.5 cm long x 1.5 cm wide x 2.0 cm high, 1 /Pkg; 10 Pkgs/Box	10
Q770407	ThinPack™ Nasal Packing, 4000Plus, Sterile 4.5 cm long x 1.0 cm wide x 2.0 cm high, 1 /Pkg; 10 Pkgs/Box	10
Q770404	Nasal Packing, 4000Plus, Sterile 8.0 cm long x 1.5 cm wide x 2.0 cm high, 1 /Pkg; 10 Pkgs/Box	10
Q770406	ThinPack™ Nasal Packing, 4000Plus, Sterile 8.0 cm long x 1.0 cm wide x 2.0 cm high, 1 /Pkg; 10 Pkgs/Box	10
Q770319	Nasal Packing, No String, 4000Plus, Sterile 8.0 cm long x 1.5 cm wide x 2.0 cm high, 1 /Pkg; 20 Pkgs/Box	20





IVALON


Nasal Packing

Nasal Packings with Airway Tube

IVALON® Nasal Packings with Airway Tube are constructed of PVAc with the 4000Plus non-stick surface.

- Integral drawstring facilitates removal

Product Code	Product Description	Per Box
Q770408	Nasal Packing with Airway Tube, 4000Plus, Sterile 4.5 cm long x 1.5 cm wide x 2.0 cm high 1 /Pkg; 10 Pkgs/Box	10
Q770405	Nasal Packing with Airway Tube, 4000Plus, Sterile 8.0 cm long x 1.5 cm wide x 2.0 cm high 1 /Pkg; 10 Pkgs/Box	10




IVALON® Epistaxis Packings

IVALON® Epistaxis Packings are constructed of PVAc.

- Small for treatment of anterior and intermediate epistaxis.
- Large for treatment of posterior and indeterminate epistaxis
- ThinPack™ for treatment of anterior epistaxis. Ideal for emergency use.
- Integral drawstring facilitates removal
- Compressed for easy insertion

Product Code	Product Description	Per Box
Q602315	Small Epistaxis Packing, Sterile 5.5 cm long x 1.5 cm wide x 2.5 cm high 1 /Pkg; 20 Pkgs/Box	20
Q602310	Large Epistaxis Packing, Sterile 10.0 cm long x 1.5 cm wide x 2.5 cm high 1 /Pkg; 10 Pkgs/Box	10
Q602330	ThinPack™ Epistaxis Packing, Sterile 10.0 cm long x 1.0 cm wide x 2.5 cm high 1 /Pkg; 10 Pkgs/Box	10



IVALON® Post-Op Sinus Packings

IVALON® Post-Op Sinus Packings are constructed of PVAc with the 4000Plus non-stick surface.

- Integral drawstring facilitates removal

Product Code	Product Description	Per Box
Q770530	Post-Op Sinus Packing, K6, 4000Plus, Sterile 3.5 cm long x 0.6 cm wide x 1.2 cm high, 2 /Pkg; 5 Pkgs/Box	10
Q770532	Post-Op Sinus Packing, K9, 4000Plus, Sterile 3.5 cm long x 0.9 cm wide x 1.2 cm high, 2 /Pkg; 5 Pkgs/Box	10

