



**Cantaert**

MedicalSupplies

T +32 9 282 24 44  
[www.cantaertmedical.be](http://www.cantaertmedical.be)


# IVALON

## Ear Wicks

### Ear Wicks

IVALON® Ear Wicks are constructed of PVAc.

Product Code	Product Description	Per Box
Q605200	Ear Wick, Sterile 9 mm diameter x 15 mm long, 1 /Pkg; 25 Pkgs/Box	25
Q605201	Ear Wick, Sterile 9 mm diameter x 15 mm long, 1 /Pkg; 50 Pkgs/Box	50
Q605202	Ear Wick, Sterile 9 mm diameter x 24 mm long, 1 /Pkg; 50 Pkgs/Box	50



### Fenestrated Ear Wicks

IVALON® Fenestrated Ear Wicks are constructed of PVAc.

Product Code	Product Description	Per Box
Q605205	Fenestrated Ear Wick, Sterile 9 mm diameter x 15 mm long, 1 /Pkg; 50 Pkgs/Box	50
Q605206	Fenestrated Ear Wick, Sterile 9 mm diameter x 24 mm long, 1 /Pkg; 50 Pkgs/Box	50



### Post-Op Ear Packing

IVALON® Post-Op Ear Packings are constructed of PVAc.

- Integral drawstring facilitates removal

Product Code	Product Description	Per Box
Q605204	Small Post-Op Ear Packing, Sterile 12 mm diameter x 24 mm long, 1 /Pkg; 10 Pkgs/Box	10
Q605203	Large Post-Op Ear Packing, Sterile 15 mm diameter x 24 mm long, 1 /Pkg; 10 Pkgs/Box	10



**Cantaert**  
MedicalSupplies

T +32 9 282 24 44  
[www.cantaertmedical.be](http://www.cantaertmedical.be)