

# Flat scale with high capacity

**seca**  
Precision for health

## seca 813

- + Solid steel frame, up to 440 lbs capacity
- + Non-slip rubber base for safe standing during measuring
- + Large LCD digits for easy reading
- + Step-off function for weighing without switching on the scale
- + More than 180 years of experience for maximum precision
- + From the global market leader for medical weighing and measuring



seca 813

## Flat scale with high capacity



### A robust, sturdy scale – with a touch of elegance

seca has designed the seca 813 with a bluetooth interface to transfer the Auto-HOLD value automatically to an external device. Thanks to its solid steel frame, the seca 813 can withstand high weight and features large LCD digits for easy reading. Furthermore, its extrawide surface and low-profile facilitate trouble-free weighing.



### Pleasant and comfortable weighing

All platforms are not created equal. For this weighing platform seca has chosen a highquality, two-component raised profile material which is soft like rubber and pleasant to the touch. But don't be fooled – this material is also tough. It fulfils two requirements at the same time: pleasant weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

Capacity	440 lbs / 200 kg
Graduation	0.2 lb / 100 g
Dimensions (WxHxD)	17 x 1.9 x 14.7 in / 433 x 47 x 373 mm
Platform (WxHxD)	15.7 x 1.8 x 11.8 in / 400 x 45 x 300 mm
Net weight	6.4 lbs / 2.9 kg
Power supply	Batteries
Functions	Auto-HOLD, Step-off, Tap-on automatic switch-on, Automatic switch-off, lbs/kg switch-over



### Paying attention to details

Weight and girth are no problem for the seca 813. You don't even have to bend over to switch it on, thanks to the unique step-off function. An issue with many scales is the display. Thanks to the 1.1 inch high LCD digits, it is extremely easy to read the measurements. You see, a scale from seca makes life more simple by paying attention to even the slightest details.

We reserve the right to make modifications.

seca north america | west  
13601 Benson Avenue  
Chino, CA 91710  
USA

seca north america | east  
1439 Stuart Engals Blvd, Suite 200  
Mount Pleasant, SC 29464  
USA

T +1 800 542 7322  
F +1 888 705 7397  
E info.us@seca.com

Precision for health  
[seca.com](http://seca.com)

5/2023