

Setpack® Toptex® lite Laparotomy Swabs X-ray detectable

Economical | safe | flexible | ready-for-use



Fields of application

- designed for immediate use in the OR
- for grasping and retaining organs and bodily tissues
- for haemostasis and for absorbing body fluids during surgical procedures

Properties

- made from absorbent cotton gauze in compliance with DIN EN 14079
- prewashed
- soft and absorbent
- 20 threads
- packaged in inner wrapping for sterile materials for fast and safe removal of the product and a swab count
- separately attached, double, self-adhesive counting cards

Product composition

100% cotton, X-ray detectable chip (polyvinyl chloride and barium sulphate)

Versions

- 4-layered and 6-layered

Setpack® Toptex® lite Laparotomy Swabs X-ray detectable

sterile

	Contents	Size (cm)	REF	Shipping unit (Pack)
white				
2	pcs., 4-layered	40 x 40	15004	14/42
5	pcs., 4-layered	40 x 40	15001	6/18
2	pcs., 6-layered	40 x 40	15012	12/36
5	pcs., 6-layered	40 x 40	15003	6/18
2	pcs., 4-layered	20 x 30	15011	28/84
5	pcs., 4-layered	20 x 30	15007	14/42
2	pcs., 4-layered	50 x 60	15014	8/24
2	pcs., 4-layered	8 x 90	15010	26/78
green				
2	pcs., 4-layered	40 x 40	15006	14/42
5	pcs., 4-layered	40 x 40	15002	6/18
2	pcs., 6-layered	40 x 40	15008	12/36
5	pcs., 6-layered	40 x 40	15000	6/18
2	pcs., 4-layered	20 x 30	15017	28/84
5	pcs., 4-layered	20 x 30	15005	14/42
2	pcs., 4-layered	50 x 60	15013	8/24
2	pcs., 6-layered	50 x 60	15018	6/18
2	pcs., 4-layered	8 x 90	15009	26/78
blue				
2	pcs., 4-layered	40 x 40	15016	14/42
5	pcs., 4-layered	40 x 40	15015	6/18

Not available in all countries