

## Actimove® TaloCast Air

### Ankle Brace

Anatomically shaped stabilizer that secures the upper ankle joint. Adjustable air cells for compression and a reliable and easy to adjust heel strap.

- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- Pre-inflated, adjustable air cells offer compression for a customized fit.
- Long-lasting, extra durable easily adjustable heel strap on the outside of the shell.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as liner and shells can be wiped with a damp cloth.
- Fits left and right foot and simplifies stock keeping.
- Available in two sizes.
- For stock keeping, long shelf life of five years.



This product is intended to stabilize the ankle joint.

Fields of application include:

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema



## Actimove® TaloCast AirGel

### Ankle Brace

Anatomically shaped stabilizer that secures the upper ankle joint. Adjustable air cells for compression and optional cold therapy for extra pain relief.

- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- Pre-inflated, adjustable air cells offer compression for a customized fit.
- Extra pain management through integrated gel pad for cold therapy.
- Long-lasting, extra durable easily adjustable heel strap on the outside of the shell.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as liner and shells can be wiped with a damp cloth.
- Fits left and right foot and simplifies stock keeping.
- Available in two sizes.
- For stock keeping, long shelf life of five years.



This product is intended to stabilize the ankle joint.

Fields of application include:

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema

Code	Size	Measurements	Unit
7311904	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
7311905	L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace



Code	Size	Measurements	Unit
7311900	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
7311901	L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace

