

# NON-SAFETY ENTERAL FEEDING DEVICES


## FUNNELED FEEDING TUBES

### PVC FEEDING TUBES

#### WITHOUT WEIGHTED TIP

The funneled PVC feeding tubes are intended to be used for short term use (up to 5 days) and feature the following characteristics:

- Closed or open tip
- Proximal end:
  - Small flexible funnel with cap with colour code for Fr 4 to Fr 10
  - Large flexible funnel without cap for Fr 12 to Fr 21
- Centimetric markings every centimetre and numerical markings every 5cm:
  - from 5cm to 35cm for codes 310 & 311
  - from 5cm to 70cm for code 312
  - from 5cm to 75cm for code 391 & 390
- Radiopaque.
- Transparent to visualize the administered solution.



Large funnel



Small funnel

Code	Tube				Funnel		Quantity
	Length cm	Fr	Tip	Lateral eyes	Type	Colour code	Box/case
310.04	40	4	Closed	2	Small	White	50/600
310.05	40	5	Closed	2	Small	Grey	50/600
310.06	40	6	Closed	2	Small	Green	50/600
310.07	40	7	Closed	2	Small	Yellow	50/600
310.08	50	8	Closed	2	Small	Blue	50/600
310.10	50	10	Closed	2	Small	Black	50/600
310.12	50	12	Closed	2	Large	-	50/600
311.04	75	4	Closed	2	Small	White	50/600
311.05	75	5	Closed	2	Small	Grey	50/600
311.06	75	6	Closed	2	Small	Green	50/600
311.08	75	8	Closed	2	Small	Blue	50/600
311.10	75	10	Closed	2	Small	Black	50/600
312.04	125	4	Closed	2	Small	White	25/600
312.06	125	6	Closed	2	Small	Green	25/600
312.07	125	7	Closed	2	Small	Yellow	25/600
312.08	125	8	Closed	2	Small	Blue	25/600
312.09	125	9	Closed	2	Small	Brown	25/600
312.10	125	10	Closed	2	Small	Black	25/600
312.12	125	12	Closed	2	Large	-	25/300
391.06	125	6	Closed	4	Small	Green	50/450
391.08	125	8	Closed	4	Small	Blue	50/450
391.10	125	10	Closed	4	Small	Black	50/300
391.12	125	12	Closed	4	Large	-	50/300
391.14	125	14	Closed	4	Large	-	50/300
391.16	125	16	Closed	4	Large	-	50/200
391.18	125	18	Closed	4	Large	-	50/200
391.21	125	21	Closed	4	Large	-	50/200
390.10	100	10	Open	4	Small	Black	50/300
390.12	100	12	Open	4	Large	-	50/300
390.14	100	14	Open	4	Large	-	50/300
390.16	100	16	Open	4	Large	-	50/300
390.18	100	18	Open	4	Large	-	50/200
390.21	100	21	Open	4	Large	-	50/200