

Electrosurgical Electrodes, Sterile, Single Use

Tungsten Microsurgical Needle Electrodes



E1650

10 units/case

Straight Microsurgical Needle

- Insulated
- Total length: 4,2 cm (1,65 in.)
- Exposed length: 2 cm (0,8 in.)
- Active length: 3,0 mm (0,120 in.)
- Tip radius: 0,06 mm (0,0023 in.)



E1652

10 units/case

Microsurgical Needle

- Insulated
- Total length: 5,4 cm (2,13 in.)
- Exposed length: 3 cm (1,2 in.)
- Active length: 3,0 mm (0,120 in.)
- Tip radius: 0,06 mm (0,0023 in.)
- 45 degree angle, 3,0 mm (0,120 in.) from tip



E1651

10 units/case

Straight Microsurgical Needle

- Insulated
- Total length: 5,4 cm (2,13 in.)
- Exposed length: 3 cm (1,2 in.)
- Active length: 3,0 mm (0,120 in.)
- Tip radius: 0,06 mm (0,0023 in.)



E1653

10 units/case

Microsurgical Needle

- Insulated
- Total length: 5,4 cm (2,13 in.)
- Exposed length: 3 cm (1,2 in.)
- Active length: 3,0 mm (0,120 in.)
- Tip radius: 0,06 mm (0,0023 in.)
- 45 degree angle, 1 cm (0,4 in.) from tip

Standard 2,4 mm (3/32 in.) shaft is designed to fit most electrosurgical pencils.

Electrosurgical Electrodes, Sterile, Single Use

Tungsten Microsurgical Needle Electrodes



E1654B

10 units/case

Microsurgical Needle With Safety Sleeve™

Not compatible with foot switching pencils or AccuVac™ smoke evacuation attachment.

- **Insulated**
- **Safety Sleeve™**
- Total length: 10,16 cm (4 in.)
- Exposed length: 9 cm (3,5 in.)
- Active length: 3,0 mm (0,120 in.)
- Tip radius: 0,06 mm (0,0023 in.)

Standard 2,4 mm (3/32 in.) shaft is designed to fit most electrosurgical pencils.



Latex
Free