


We have a range of lightweight versions of the most commonly used anaesthetic breathing systems in the Mapleson classifications.

### Mapleson A and C breathing systems




system illustrated

Code	Description	Box Qty.
2105000	Adult, Mapleson A Magill anaesthetic breathing system with 2L bag, ≥1.6m	15
2106000	Adult, Mapleson A Parallel Lack breathing system with 2L bag, ≥1.6m	10
2108000 (S*)	Adult, Mapleson C bagging system with 2L bag, ≥1.8m	15 (20*)

Make an enquiry

### Mapleson D and Deluxe Bain coaxial breathing systems



system illustrated

Code	Description	Box Qty.
2110000	Mapleson D Bain coaxial breathing system for use with Bain adaptors (no valve or bag), ≥1.6m	20
2115000	Mapleson D deluxe Bain coaxial breathing system with 2L bag, ≥1.6m	30
2115100	Silver Knight™ anti-microbial Mapleson D deluxe Bain coaxial breathing system with 2L bag, ≥1.6m	30
2116000	Mapleson D deluxe Bain coaxial breathing system with 2L bag, ≥2.8m	10
2117000	Mapleson D deluxe Bain coaxial breathing system with 2L bag, ≥5.4m	5

Make an enquiry

Silver Knight

Sterile option available

(S\*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S (sterile box quantity is shown in brackets).



Contact us