

# Telacomp®



## Description

The Telacomp® gauze swabs have cut edges folded in and have an interwoven Telatrast® X-ray detectable thread. The sterile Telacomp® gauze swabs are supplied in a double peel pack with double control cards. Also available as non-sterile version.

## Uses

For staunching bleeding and wound treatment as well as for gripping and retaining organs and body tissue.



Product	Size	Package contents	Article number
<b>Telacomp® sterile</b> Double sterile packaging 12-ply, 17 threads	7.5 x 7.5 cm	12 x 10 Pieces	452 133
	10 x 10 cm	16 x 10 Pieces	452 135
	10 x 10 cm	12 x 20 Pieces	452 235
	10 x 20 cm	12 x 10 Pieces	452 137
<b>Telacomp® steriles</b> Double sterile packaging 16-ply, 17 threads	7.5 x 7.5 cm	10 x 10 Pieces	452 143
	10 x 10 cm	14 x 10 Pieces	452 145
<b>Telacomp® non-sterile</b> 8-ply, 17 threads	10 x 10 cm	100 Pieces	401 925
<b>Telacomp® non-sterile</b> 12-ply, 17 threads	7.5 x 7.5 cm	100 Pieces	401 933
	10 x 10 cm	100 Pieces	401 935
	10 x 20 cm	100 Pieces	401 937
<b>Telacomp® non-sterile</b> 16-ply, 17 threads	7.5 x 7.5 cm	100 Pieces	401 943
	10 x 10 cm	100 Pieces	401 945