


HOSPITAL CARE

Terumo Syringe with Needle

3-part Syringe with Needle Product Information



Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

- The integral gasket ensures a firm and distinct end-stop
- Clear graduated scale allows accurate dosage
- Transparent barrel for visualisation
- Finger grips to support easy handling
- Lubricated needle and precisely shaped needle point

Description		Raw Materials
Syringe		
1	Barrel	Polypropylene
2	Plunger	Polypropylene
3	Gasket	Thermoplastic elastomer (for ,SS'-codes) Polyisoprene rubber (for ,MDSS'-codes)
	Lubricant	Silicone oil
Needle		
4	Hub	Polypropylene + master batch
5	Cannula	Stainless steel
	Adhesive	Epoxy glue
	Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.

Sterilization

With radiation (e-beam) (for ,SS'-codes) or ethylene oxide (for ,MDSS'-codes) to a Sterility Assurance Level of at least 10^{-6}

Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper (for ,MDSS'-codes) or embossed in the film of the blister pack (for ,SS'-codes)
- Do not store at extreme temperature and humidity

Pre-connected with
a hypodermic needle

Product Information Terumo Syringe with Needle



Product Range

Product Code	Nominal Capacity	Cannula ø		Cannula Length		Replacement for
MDSS01S2516E	1 ml	25G	0.5 mm	5/8 inch	16 mm	SS+01T25161
MDSS01S2613E	1 ml	26G	0.45 mm	1/2 inch	13 mm	SS+01T26131
MDSS01S2713E	1 ml	27G	0.4 mm	1/2 inch	13 mm	–
SS*02SE21251	2.5 ml	21G	0.8 mm	1 inch	25 mm	–
SS*02SE22381	2.5 ml	22G	0.7 mm	1 1/2 inch	38 mm	–
SS*02SE23321	2.5 ml	23G	0.6 mm	1 1/4 inch	32 mm	–
MDSS03S2116E	3 ml	21G	0.8 mm	5/8 inch	16 mm	SS*02SE21161
MDSS03S2125E	3 ml	21G	0.8 mm	1 inch	25 mm	** SS*02SE21251
MDSS03S2138E	3 ml	21G	0.8 mm	1 1/2 inch	38 mm	SS*02SE21381
MDSS03S2325E	3 ml	23G	0.6 mm	1 inch	25 mm	SS*02SE23251
MDSS03S2516E	3 ml	25G	0.5 mm	5/8 inch	16 mm	SS*02SE25161
MDSS05S2138E	5 ml	21G	0.8 mm	1 1/2 inch	38 mm	SS+05S21381
MDSS05S2238E	5 ml	22G	0.7 mm	1 1/2 inch	38 mm	SS+05S22381
MDSS05S2332E	5 ml	23G	0.6 mm	1 1/4 inch	32 mm	SS+05S23321
MDSS10S2038E	10 ml	20G	0.9 mm	1 1/2 inch	38 mm	SS+10S20381
MDSS10S2138E	10 ml	21G	0.8 mm	1 1/2 inch	38 mm	SS+10S21381
MDSS10S2238E	10 ml	22G	0.7 mm	1 1/2 inch	38 mm	SS+10S22381

CE
0197



STERILE EO
for ,MDSS'-codes

STERILE R
for ,SS'-codes

** as from January 2023



Packaging

Unit Box	Shipping Carton		
	Quantity	Gross Weight	Dimensions (L x W x H)
100	2400	12.50 kg	547 x 425 x 405 mm
100	2400	12.50 kg	547 x 425 x 405 mm
100	2400	12.50 kg	547 x 425 x 405 mm
100	1800	9.70 kg	522 x 389 x 439 mm
100	1800	9.70 kg	522 x 389 x 439 mm
100	1800	9.70 kg	522 x 389 x 439 mm
100	2400	16.80 kg	580 x 425 x 485 mm
100	2400	16.80 kg	580 x 425 x 485 mm
100	2400	16.80 kg	580 x 425 x 485 mm
100	2400	16.80 kg	580 x 425 x 485 mm
100	2400	16.80 kg	580 x 425 x 485 mm
100	1800	16.60 kg	575 x 430 x 490 mm
100	1800	16.60 kg	575 x 430 x 490 mm
100	1800	16.60 kg	575 x 430 x 490 mm
100	1200	14.60 kg	625 x 458 x 480 mm
100	1200	14.60 kg	625 x 458 x 480 mm
100	1200	14.60 kg	625 x 458 x 480 mm

Single unit: Blister pack
(film/film for ,SS'-codes and film/paper for ,MDSS'-codes)



Regulatory Status

CE Certification –
Class IIa medical device

Notified body
TÜV Rheinland 0197

Manufacturer
for ,SS'-codes TERUMO CORPORATION,
2-44-1 Hatagaya, Shibuya-Ku,
Toyko 151-0072, Japan
for ,MDSS'-codes Yangzhou Medline
Industry Co., Ltd. (China)

EC representative
for ,SS'-codes TERUMO EUROPE N.V. –
Leuven, Belgium
for ,MDSS'-code Shanghai International
Corp. GmbH (Europe) – Hamburg,
Germany

Distributor
TERUMO EUROPE N.V.,
Interleuvenlaan 40, 3001 Leuven, Belgium

Precautions

- Do not use if package is damaged
- Do not re-use
- Dispose of sharps after single use according to facility protocol to avoid risk of infection

Always consult the product label for a complete overview of warnings, cautions and/or precautions prior to actual use.

Choose from a range of volumes from 1 ml to 10 ml

All syringe sizes are available with preconnected needles in different gauges.

Your Terumo partner



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

Please note that some products may not be available in all regions, check with your local Terumo contact.