

# Freezpen is a powerful, precise cryosurgery instrument for the removal of benign skin lesions.

#### Simply powerful

High power (50-55 bar) microjet achieves in-tissue temperature of -89°C/-128°F, for complete tissue necrosis in seconds.

# Simply precise

Ultrafine sharp jet for total precision and no damage to healthy tissue. Range of applicators.

#### Simple procedure

Assess the lesion, choose the right applicator, point, freeze, wipe, freeze again. Done!

# Simply painless

No pain, no cuts or bleeding, no anesthesia - even for children. Slight sting signals when treatment is complete.

## Simply perfect

Perfect results with no scarring. Healthy tissue is unaffected; frozen tissue disappears in a few weeks.

## Simply practical

Easy handling and transport. Handy disposable cartridges; store at room temperature.

# It's simply good business!

Each treatment takes just seconds. Return on investment in 10 weeks at just six treatments/week.



# Freezpen can remove most common skin pathologies in just a few seconds.





**Freezpen 8g** instrument set - 2 applicators: Blue and white applicators. User manual and tweezers for filter placement.



Box with 21 N<sub>2</sub>O cartridges of 8g each carton with 12 boxes.



**Freezpen 16g** instrument set - 3 applicators: Blue, white and green applicators. User manual and tweezers for filter placement.



Box with 12 N<sub>2</sub>O cartridges of 16g each carton with 12 boxes.



**Freezpen 35g** instrument set - 3 applicators: Green, orange long and white long applicators. User manual and tweezers for filter placement.



Box with 6 N<sub>2</sub>O cartridges of 35g each carton with 12 boxes.



