

BD Insyte™ Autoguard™ peripheral safety straight IV catheter

- Push-button shielding technology proven to reduce needlestick injuries by 95 % vs. a non-safety catheter²
- Provides safe handling of a withdrawn needle without reducing its ease of insertion⁶
- Pressing the push-button releases the spring so that the needle/flashback chamber retracts into the safety barrel
- BD Luer-Lok™ compatible
- Includes all the benefits of BD Vialon™ (see page 7)
- BD Instaflash™ needle technology immediately confirms vessel entry at the point of insertion and is designed to improve insertion success
- Sterile, single use



Straight

Gauge	Catheter Length (mm)	Outer diameter (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference
24 G	19	0.7	20	Yellow	50 / 200	381812
22 G	25	0.9	35	Blue	50 / 200	381823
20 G	25	1.1	65	Pink	50 / 200	381833
20 G	30	1.1	60	Pink	50 / 200	381834
20 G	48	1.1	55	Pink	50 / 200	381837
18 G	30	1.3	105	Green	50 / 200	381844
18 G	48	1.3	95	Green	50 / 200	381847
16 G	30	1.7	220	Grey	50 / 200	381854
16 G	45	1.7	205	Grey	50 / 200	381857
14 G	45	2.1	330	Orange	50 / 200	381867

Winged

Gauge	Catheter Length (mm)	Outer diameter (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference
24 G	19	0.7	20	Yellow	50 / 200	381912
22 G	25	0.9	35	Blue	50 / 200	381923
20 G	30	1.1	65	Pink	50 / 200	381933
20 G	30	1.1	60	Pink	50 / 200	381934
20 G	48	1.1	55	Pink	50 / 200	381937
18 G	30	1.3	105	Green	50 / 200	381944
18 G	48	1.3	95	Green	50 / 200	381947
16 G	30	1.7	220	Grey	50 / 200	381954
16 G	45	1.7	205	Grey	50 / 200	381957