

BD™ Blunt Fill and Blunt Filter Needles

- Fill needle with blunt tip
- 45° bevel avoids coring and designed to minimise the risk of needlestick injury
- Red hub and shield for easy identification as fill needle
- Also available with 5 micron filter for effective filtration of glass, metal, rubber and other contaminant particles
- Sterile, single use



BD™ Blunt Fill Needles

Product Code	Gauge	Length (mm)	Length (inch)	Colour code	Units/ Box	NHS Supply Chain Code
303129†	18 G	40	1 1/2"	Red	100	FTR1988
305181	18 G	25	1"	Red	100	-

† Product also available in bulk non sterile.

BD™ Blunt Fill Needle with 5 Micron Filter

Product Code	Gauge	Length (mm)	Length (inch)	Colour code	Units/ Box	NHS Supply Chain Code
305211†	18 G	40	1 1/2"	Purple	100	FTR436

† Product also available in bulk non sterile.

BD Emerald™ with BD™ Blunt Fill Needle

- BD Emerald™ with BD™ Blunt Fill Needle pre-attached for safe and efficient medication preparation process
- Sterile, single use



BD Emerald™ with BD™ Blunt Fill Needle

Product Code	Description	Gauge	Capacity (ml)	Length (mm)	Length (inch)	Units/ Box	NHS Supply Chain Code
303221	BD Emerald™ with BD™ Blunt Fill Needle	18 G	2	40	1 1/2"	100	-
303139	BD Emerald™ with BD™ Blunt Fill Needle	18 G	5	40	1 1/2"	100	-
303140	BD Emerald™ with BD™ Blunt Fill Needle	18 G	10	40	1 1/2"	100	-